

BOA TOURNIQUET

HYPER CARRY®

UPPER AND LOWER EXTREMITY TOURNIQUET

Hemostatic Tourniquet

In situations requiring emergency intervention, it is vital to control bleeding immediately. The hemostatic tourniquet BOA which is developed under the brand name family as Hyper Carry, is an effective, durable and reliable solution designed for military personnel, healthcare teams, security forces and first aid officers.

**TURK
[PATENT]**
TÜRK PATENT VE MARKA KURUMU
DESIGN REGISTRATION
2025/031970



BOA TOURNIQUET

- Fast One-Handed Application
- Strong and Durable Material
- Safe Locking Mechanism
- Compact and Lightweight Design
- Compliance with Military and Medical Standards

for further information
hypercarry.com.tr

FAST AND EFFECTIVE INTERVENTION

Rapid control of severe bleeding is vital on the battlefield. Extremity tourniquets are life-saving tools used to manage major arterial injuries, particularly in the arms and legs. In military first aid interventions, a properly placed and tightened tourniquet helps to minimize blood loss and keeps the casualty stable until evacuation. When applied correctly and promptly, a tourniquet plays a critical role in reducing preventable deaths.

① Time Strap
It is designed to write the time of the tourniquet application.

③ BOA Windlass
Made from durable polyamide material.

⑤ BOA Plate
Flexible and durable structure, able to take the shape of the application surface.

⑦ BOA Buckle
It allows the webbing to wrap around the extremity and helps tighten it.

② Safety Orange Tip
Coloured tip for easy application.

④ C Lock
Slot designed for fixing the windlass.

⑥ Webbing
Webbing made of durable material with special velcro sewn on that can stick to itself. Can be applied on wet and dirty surfaces.

PRODUCT SPECIFICATIONS

Packaged Size & Weight

Length: 2.5 in (6.5 cm)
Width: 2 in (5 cm)
Height: 6.3 in (16 cm)
Weight: 3.175 oz (90 g)

Models

Ref No: TQ.0001 (Black)



Hyper Carry
Dalğıç Savunma Sanayi ve Tic. Ltd. Şti.



Fatih Mah. Emek Cd. No: 27/A
Pursaklar Ankara Türkiye
Made in Türkiye



Mail: info@hypercarry.com.tr
Phone: +90 312 874 50 10
Fax: +90 312 874 50 11



for further information
hypercarry.com.tr